

Product Datasheet

Porcine Interleukin 15 ELISA Kit (orb1211961)

Description	Porcine Interleukin 15, IL-15 ELISA Kit
Reactivity	Porcine
Range	0.3-90ng/L
Tested Applications	ELISA
Target	IL15
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,437
Uniprot ID	Q95253
Solubility (25°C)	96ng/L
Sensitivity	0.152ng/L
Entrez	397683
Expiration Date	Please enquire.