

Product Datasheet

Porcine Interleukin 6 ELISA Kit (orb1211963)

Description	Porcine Interleukin 6, IL-6 ELISA Kit
Reactivity	Porcine
Range	20-6000ng/L
Tested Applications	ELISA
Target	IL6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,953
Uniprot ID	P26893
Solubility (25°C)	6400ng/L
Sensitivity	10.06ng/L
Entrez	100628202399500
Expiration Date	Please enquire.