

Product Datasheet

Mouse Glutamate Receptor Ionotropic Nmdar2B ELISA Kit (orb1211996)

Description	Mouse Glutamate Receptor Ionotropic Nmdar2B, GRIN2B ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	GRIN2B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	165,959
Uniprot ID	Q01097
Solubility (25°C)	32ng/ml
Sensitivity	0.032ng/ml
Entrez	14812
Expiration Date	Please enquire.