

Product Datasheet

Mouse Atpase H⁺/K⁺ Transporting Subunit Alpha ELISA Kit (orb1212007)

Description	Mouse Atpase H ⁺ /K ⁺ Transporting Subunit Alpha, ATP4A ELISA Kit
Reactivity	Mouse
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	SLC51A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,759
Uniprot ID	Q8R000
Solubility (25°C)	12.8ng/ml
Sensitivity	0.015ng/ml
Entrez	106407
Expiration Date	Please enquire.