

Product Datasheet

Mouse Triggering Receptor Expressed On Myeloid Cells 2 ELISA Kit (orb1212033)

Description	Mouse Triggering Receptor Expressed On Myeloid Cells 2, TREM2 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	TREM2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,527
Uniprot ID	Q99NH8
Solubility (25°C)	24ng/mL
Sensitivity	0.024ng/mL
Entrez	83433
Expiration Date	Please enquire.