

Product Datasheet

Mouse Syndecan 4 ELISA Kit (orb1212043)

Description	Mouse Syndecan 4, SDC4 ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	SDC4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,482
Uniprot ID	035988
Solubility (25°C)	480ng/L
Sensitivity	0.56ng/L
Entrez	20971
Expiration Date	Please enquire.