

Product Datasheet

Mouse Spondin-1 ELISA Kit (orb1212046)

Description	Mouse Spondin-1, SPON1 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	SPON1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	90,821
Uniprot ID	Q8VCC9
Solubility (25°C)	24ng/mL
Sensitivity	0.023ng/mL
Entrez	233744
Expiration Date	Please enquire.