

Product Datasheet

Mouse Selenoprotein P ELISA Kit (orb1212053)

Description	Mouse Selenoprotein P, SEPP1 ELISA Kit
Reactivity	Mouse
Range	0.1-40ug/mL
Tested Applications	ELISA
Target	SELENOP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,706
Uniprot ID	P70274
Solubility (25°C)	48µg/mL
Sensitivity	0.056µg/mL
Entrez	20363
Expiration Date	Please enquire.