

Product Datasheet

Mouse Secretoglobin Family 3A Member 2 ELISA Kit (orb1212056)

Description	Mouse Secretoglobin Family 3A Member 2, SCGB3A2 ELISA Kit
Reactivity	Mouse
Range	5-1000ng/L
Tested Applications	ELISA
Target	SCGB3A2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	9,819
Uniprot ID	Q920H1
Solubility (25°C)	1280ng/L
Sensitivity	2.15ng/L
Entrez	117158
Expiration Date	Please enquire.