

Product Datasheet

Mouse R-spondin-2 ELISA Kit (orb1212058)

Description	Mouse R-spondin-2, RSPO2 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	RSPO2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,276
Uniprot ID	Q8BFU0
Solubility (25°C)	4800ng/L
Sensitivity	9.95ng/L
Entrez	239405
Expiration Date	Please enquire.