

Product Datasheet

Mouse Reelin ELISA Kit (orb1212064)

Description	Mouse Reelin, RELN ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	RELN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	387,495
Uniprot ID	Q60841
Solubility (25°C)	960ng/L
Sensitivity	1.54ng/L
Entrez	19699
Expiration Date	Please enquire.