

Product Datasheet

Mouse Protein S100-A4 ELISA Kit (orb1212068)

Description	Mouse Protein S100-A4, S100A4 ELISA Kit
Reactivity	Mouse
Range	30-9000ng/L
Tested Applications	ELISA
Target	S100A4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,721
Uniprot ID	P07091
Solubility (25°C)	9600ng/L
Sensitivity	14.16ng/L
Entrez	20198
Expiration Date	Please enquire.