

## Product Datasheet

### Mouse Protein Argonaute-2 ELISA Kit (orb1212073)

|                            |                                                |
|----------------------------|------------------------------------------------|
| <b>Description</b>         | Mouse Protein Argonaute-2,<br>EIF2C2 ELISA Kit |
| <b>Reactivity</b>          | Mouse                                          |
| <b>Range</b>               | 20-4500ng/L                                    |
| <b>Tested Applications</b> | ELISA                                          |
| <b>Target</b>              | AGO2                                           |
| <b>Storage</b>             | 2-8°C                                          |
| <b>Note</b>                | For research use only                          |
| <b>Application notes</b>   | Sandwich                                       |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates         |
| <b>MW</b>                  | 97,304                                         |
| <b>Uniprot ID</b>          | <b>Q8CJG0</b>                                  |
| <b>Solubility (25°C)</b>   | 4800ng/L                                       |
| <b>Sensitivity</b>         | 9.22ng/L                                       |
| <b>Entrez</b>              | <b>239528</b>                                  |
| <b>Expiration Date</b>     | Please enquire.                                |