

Product Datasheet

Mouse Prokineticin-2 ELISA Kit (orb1212075)

Description	Mouse Prokineticin-2, PROK2 ELISA Kit
Reactivity	Mouse
Range	7-1500ng/L
Tested Applications	ELISA
Target	PROK2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,185
Uniprot ID	Q9QXU7
Solubility (25°C)	1600ng/L
Sensitivity	3.43ng/L
Entrez	50501
Expiration Date	Please enquire.