

Product Datasheet

Mouse Plastin-2 ELISA Kit (orb1212077)

Description	Mouse Plastin-2, LCP1 ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	LCP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	70,149
Uniprot ID	Q61233
Solubility (25°C)	3200ng/L
Sensitivity	7.36ng/L
Entrez	18826
Expiration Date	Please enquire.