

Product Datasheet

Mouse Osteoclast-associated Immunoglobulin-like Receptor ELISA Kit (orb1212094)

Description	Mouse Osteoclast-associated Immunoglobulin-like Receptor, OSCAR ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	OSCAR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,726
Uniprot ID	Q8VBT3
Solubility (25°C)	4800ng/L
Sensitivity	9.38ng/L
Entrez	232790
Expiration Date	Please enquire.