

Product Datasheet

Mouse N-Myc-Interactor ELISA Kit (orb1212098)

Description	Mouse N-Myc-Interactor, NMI ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	NMI
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,236
Uniprot ID	035309
Solubility (25°C)	2400ng/L
Sensitivity	5.19ng/L
Entrez	64685
Expiration Date	Please enquire.