

Product Datasheet

Mouse Neurturin ELISA Kit (orb1212102)

Description	Mouse Neurturin, NRTN ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	NRTN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,219
Uniprot ID	P97463
Solubility (25°C)	4800ng/L
Sensitivity	9.87ng/L
Entrez	18188
Expiration Date	Please enquire.