

Product Datasheet

Mouse Neuregulin-4 ELISA Kit (orb1212109)

Description	Mouse Neuregulin-4, NRG4 ELISA Kit
Reactivity	Mouse
Range	1-400ng/mL
Tested Applications	ELISA
Target	NRG4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,744
Uniprot ID	Q9WTX4
Solubility (25°C)	480ng/mL
Sensitivity	0.54ng/mL
Entrez	83961
Expiration Date	Please enquire.