

Product Datasheet

Mouse Myc Proto-oncogene Protein ELISA Kit (orb1212116)

Description	Mouse Myc Proto-oncogene Protein, MYC ELISA Kit
Reactivity	Mouse
Range	20-3500ng/L
Tested Applications	ELISA
Target	MYC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	48,971
Uniprot ID	P01108
Solubility (25°C)	4000ng/L
Sensitivity	9.05ng/L
Entrez	17869
Expiration Date	Please enquire.