

Product Datasheet

Mouse Mitochondrial Brown Fat Uncoupling Protein 1 ELISA Kit (orb1212122)

Description	Mouse Mitochondrial Brown Fat Uncoupling Protein 1, UCP1 ELISA Kit
Reactivity	Mouse
Range	20-6000ng/L
Tested Applications	ELISA
Target	UCP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,248
Uniprot ID	P12242
Solubility (25°C)	6400ng/L
Sensitivity	9.27ng/L
Entrez	22227
Expiration Date	Please enquire.