

Product Datasheet

Mouse Macrophage Scavenger Receptor Types 1 and 2 ELISA Kit (orb1212139)

Description	Mouse Macrophage Scavenger Receptor Types 1 and 2, MSR1 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	MSR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,170
Uniprot ID	P30204
Solubility (25°C)	24ng/mL
Sensitivity	0.024ng/mL
Entrez	20288
Expiration Date	Please enquire.