

## Product Datasheet

# Mouse Macrophage Mannose Receptor 1 ELISA Kit (orb1212140)

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Description</b>         | Mouse Macrophage Mannose Receptor 1, MRC1 ELISA Kit |
| <b>Reactivity</b>          | Mouse                                               |
| <b>Range</b>               | 0.3-90ng/mL                                         |
| <b>Tested Applications</b> | ELISA                                               |
| <b>Target</b>              | MRC1                                                |
| <b>Storage</b>             | 2-8°C                                               |
| <b>Note</b>                | For research use only                               |
| <b>Application notes</b>   | Sandwich                                            |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates              |
| <b>MW</b>                  | 164,981                                             |
| <b>Uniprot ID</b>          | <b>Q61830</b>                                       |
| <b>Solubility (25°C)</b>   | 96ng/mL                                             |
| <b>Sensitivity</b>         | 0.152ng/mL                                          |
| <b>Entrez</b>              | <b>17533</b>                                        |
| <b>Expiration Date</b>     | Please enquire.                                     |