

Product Datasheet

Mouse Lysyloxidase Homolog 2 ELISA Kit (orb1212141)

Description	Mouse Lysyloxidase Homolog 2, LOXL2 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	LOXL2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	87,003
Uniprot ID	P58022
Solubility (25°C)	2400ng/L
Sensitivity	4.55ng/L
Entrez	94352
Expiration Date	Please enquire.