

Product Datasheet

Mouse Keratin ELISA Kit (orb1212156)

Description	Mouse Keratin, Type 1 Cytoskeletal 18, KRT18 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	KRT18
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,538
Uniprot ID	P05784
Solubility (25°C)	24ng/mL
Sensitivity	0.019ng/mL
Entrez	16668
Expiration Date	Please enquire.