

Product Datasheet

Mouse Kelch-like Ech-associated Protein 1 ELISA Kit (orb1212157)

Description	Mouse Kelch-like Ech-associated Protein 1, KEAP1 ELISA Kit
Reactivity	Mouse
Range	2-600ng/L
Tested Applications	ELISA
Target	KEAP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	69,553
Uniprot ID	Q9Z2X8
Solubility (25°C)	640ng/L
Sensitivity	0.97ng/L
Entrez	50868
Expiration Date	Please enquire.