

Product Datasheet

Mouse Interleukin-36 Beta ELISA Kit (orb1212160)

Description	Mouse Interleukin-36 Beta, IL1F8 ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	IL36B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,878
Uniprot ID	Q9D6Z6
Solubility (25°C)	960ng/L
Sensitivity	1.47ng/L
Entrez	69677
Expiration Date	Please enquire.