

Product Datasheet

Mouse Interleukin-36 Alpha ELISA Kit (orb1212161)

Description	Mouse Interleukin-36 Alpha, IL1F6 ELISA Kit
Reactivity	Mouse
Range	6-1800ng/L
Tested Applications	ELISA
Target	IL36A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,015
Uniprot ID	Q9JLA2
Solubility (25°C)	1920ng/L
Sensitivity	3.13ng/L
Entrez	54448
Expiration Date	Please enquire.