

Product Datasheet

Mouse Interleukin-27 Subunit Beta ELISA Kit (orb1212163)

Description	Mouse Interleukin-27 Subunit Beta, EBI3 ELISA Kit
Reactivity	Mouse
Range	0.5-200ng/L
Tested Applications	ELISA
Target	EBI3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,353
Uniprot ID	035228
Solubility (25°C)	240ng/L
Sensitivity	0.26ng/L
Entrez	50498
Expiration Date	Please enquire.