

Product Datasheet

Mouse Interleukin-24 ELISA Kit (orb1212165)

Description	Mouse Interleukin-24, IL24 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	IL24
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,812
Uniprot ID	Q925S4
Solubility (25°C)	4800ng/L
Sensitivity	9.27ng/L
Entrez	93672
Expiration Date	Please enquire.