

Product Datasheet

Mouse Interleukin-22 Receptor Subunit Alpha-2 ELISA Kit (orb1212168)

Description	Mouse Interleukin-22 Receptor Subunit Alpha-2, IL22RA2 ELISA Kit
Reactivity	Mouse
Range	20-6000ng/L
Tested Applications	ELISA
Target	IL22RA2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,613
Uniprot ID	Q80XF5
Solubility (25°C)	6400ng/L
Sensitivity	9.15ng/L
Entrez	237310
Expiration Date	Please enquire.