

Product Datasheet

Mouse Interferon Regulatory Factor 4 ELISA Kit (orb1212181)

Description	Mouse Interferon Regulatory Factor 4, IRF4 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	IRF4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,577
Uniprot ID	Q64287
Solubility (25°C)	2400ng/L
Sensitivity	5.34ng/L
Entrez	16364
Expiration Date	Please enquire.