

Product Datasheet

Mouse Interferon Regulatory Factor 1 ELISA Kit (orb1212182)

Description	Mouse Interferon Regulatory Factor 1, IRF1 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	IRF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,319
Uniprot ID	P15314
Solubility (25°C)	4800ng/L
Sensitivity	9.51ng/L
Entrez	16362
Expiration Date	Please enquire.