

Product Datasheet

Mouse Inter-alpha-trypsin Inhibitor Heavy Chain H3 ELISA Kit (orb1212184)

Description	Mouse Inter-alpha-trypsin Inhibitor Heavy Chain H3, ITIH3 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/mL
Tested Applications	ELISA
Target	ITIH3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	99,358
Uniprot ID	Q61704
Solubility (25°C)	2400ng/mL
Sensitivity	5.44ng/mL
Entrez	16426
Expiration Date	Please enquire.