

Product Datasheet

Mouse Integrin Beta-3 ELISA Kit (orb1212185)

Description	Mouse Integrin Beta-3, ITGB3 ELISA Kit
Reactivity	Mouse
Range	30-9000ng/L
Tested Applications	ELISA
Target	ITGB3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	86,738
Uniprot ID	054890
Solubility (25°C)	9600ng/L
Sensitivity	14.09ng/L
Entrez	16416
Expiration Date	Please enquire.