

Product Datasheet

Mouse Integrin Beta-1 ELISA Kit (orb1212186)

Description	Mouse Integrin Beta-1, ITGB1 ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	ITGB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	88,231
Uniprot ID	P09055
Solubility (25°C)	960ng/L
Sensitivity	1.65ng/L
Entrez	16412
Expiration Date	Please enquire.