

Product Datasheet

Mouse Integrin Alpha-2B ELISA Kit (orb1212187)

Description	Mouse Integrin Alpha-2B, ITGA2B ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	ITGA2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	128,955
Uniprot ID	Q62469
Solubility (25°C)	24ng/mL
Sensitivity	0.024ng/mL
Entrez	16398
Expiration Date	Please enquire.