

Product Datasheet

Mouse Insulin-like Peptide InsI5 ELISA Kit (orb1212188)

Description	Mouse Insulin-like Peptide InsI5, INSL5 ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	INSL5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,524
Uniprot ID	Q9WUG6
Solubility (25°C)	480ng/L
Sensitivity	0.57ng/L
Entrez	23919
Expiration Date	Please enquire.