

Product Datasheet

Mouse Inducible T-cell Costimulator ELISA Kit (orb1212189)

Description	Mouse Inducible T-cell Costimulator, ICOS ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	ICOS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,709
Uniprot ID	Q9WVS0
Solubility (25°C)	24ng/mL
Sensitivity	0.024ng/mL
Entrez	54167
Expiration Date	Please enquire.