

Product Datasheet

Mouse Indoleamine 2 ELISA Kit (orb1212190)

Description	Mouse Indoleamine 2, 3-dioxygenase 1, IDO1 ELISA Kit
Reactivity	Mouse
Range	30-9000ng/L
Tested Applications	ELISA
Target	IDO1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	45,641
Uniprot ID	P28776
Solubility (25°C)	9600ng/L
Sensitivity	14.32ng/L
Entrez	15930
Expiration Date	Please enquire.