

Product Datasheet

Mouse Hydroxyacidoxidase 1 ELISA Kit (orb1212191)

Description	Mouse Hydroxyacidoxidase 1, HAO1 ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	HAO1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,001
Uniprot ID	Q9WU19
Solubility (25°C)	960ng/L
Sensitivity	1.52ng/L
Entrez	15112
Expiration Date	Please enquire.