

Product Datasheet

Mouse Homer Protein Homolog 1 ELISA Kit (orb1212194)

Description	Mouse Homer Protein Homolog 1, HOMER1 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	HOMER1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,413
Uniprot ID	Q9Z2Y3
Solubility (25°C)	4800ng/L
Sensitivity	9.79ng/L
Entrez	26556
Expiration Date	Please enquire.