

Product Datasheet

Mouse Histone Deacetylase 6 ELISA Kit (orb1212196)

Description	Mouse Histone Deacetylase 6, HDAC6 ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	HDAC6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	125,787
Uniprot ID	Q9Z2V5
Solubility (25°C)	3200ng/L
Sensitivity	7.34ng/L
Entrez	15185
Expiration Date	Please enquire.