


Product Datasheet

Mouse High Mobility Group Protein B2 ELISA Kit (orb1212197)

Description	Mouse High Mobility Group Protein B2, HMGB2 ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	HMGB2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,162
Uniprot ID	P30681
Solubility (25°C)	48ng/mL
Sensitivity	0.056ng/mL
Entrez	97165
Expiration Date	Please enquire.



Standard Curve of Mouse High Mobility Group Protein B2 ELISA Kit