

Product Datasheet

Mouse Heterogeneous Nuclear Ribonucleoprotein A3 ELISA Kit (orb1212198)

Description	Mouse Heterogeneous Nuclear Ribonucleoprotein A3, HNRNPA3 ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	HNRNPA3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,652
Uniprot ID	Q8BG05
Solubility (25°C)	3200ng/L
Sensitivity	7.63ng/L
Entrez	229279
Expiration Date	Please enquire.