

Product Datasheet

Mouse Hepatocyte Growth Factor Activator ELISA Kit (orb1212199)

Description	Mouse Hepatocyte Growth Factor Activator, HGFAC ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	HGFAC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	70,568
Uniprot ID	Q9R098
Solubility (25°C)	4800ng/L
Sensitivity	9.78ng/L
Entrez	54426
Expiration Date	Please enquire.