

## Product Datasheet

### Mouse Heparan Sulfate Glucosamine 3-O-sulfotransferase 2 ELISA Kit (orb1212200)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Heparan Sulfate<br>Glucosamine 3-O-sulfotransferase 2, HS3ST2<br>ELISA Kit |
| <b>Reactivity</b>          | Mouse                                                                            |
| <b>Range</b>               | 20-4500ng/L                                                                      |
| <b>Tested Applications</b> | ELISA                                                                            |
| <b>Target</b>              | HS3ST2                                                                           |
| <b>Storage</b>             | 2-8°C                                                                            |
| <b>Note</b>                | For research use only                                                            |
| <b>Application notes</b>   | Sandwich                                                                         |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                                           |
| <b>MW</b>                  | 41,378                                                                           |
| <b>Uniprot ID</b>          | <a href="#">Q673U1</a>                                                           |
| <b>Solubility (25°C)</b>   | 4800ng/L                                                                         |
| <b>Sensitivity</b>         | 9.31ng/L                                                                         |
| <b>Entrez</b>              | <a href="#">195646</a>                                                           |
| <b>Expiration Date</b>     | Please enquire.                                                                  |