

## Product Datasheet

### Mouse Heat-Stable Enterotoxin Receptor ELISA Kit (orb1212202)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | Mouse Heat-Stable Enterotoxin Receptor, GUCY2C ELISA Kit |
| <b>Reactivity</b>          | Mouse                                                    |
| <b>Range</b>               | 20-4500ng/L                                              |
| <b>Tested Applications</b> | ELISA                                                    |
| <b>Target</b>              | GUCY2C                                                   |
| <b>Storage</b>             | 2-8°C                                                    |
| <b>Note</b>                | For research use only                                    |
| <b>Application notes</b>   | Sandwich                                                 |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                   |
| <b>MW</b>                  | 123,232                                                  |
| <b>Uniprot ID</b>          | <b>Q3UWA6</b>                                            |
| <b>Solubility (25°C)</b>   | 4800ng/L                                                 |
| <b>Sensitivity</b>         | 9.67ng/L                                                 |
| <b>Entrez</b>              | <b>14917</b>                                             |
| <b>Expiration Date</b>     | Please enquire.                                          |