

Product Datasheet

Mouse Heat Shock Protein Hsp 90-beta ELISA Kit (orb1212203)

Description	Mouse Heat Shock Protein Hsp 90-beta, HSP90AB1 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	HSP90AB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	83,281
Uniprot ID	P11499
Solubility (25°C)	2400ng/L
Sensitivity	5.28ng/L
Entrez	15516
Expiration Date	Please enquire.