

Product Datasheet

Mouse Glutathione Peroxidase 3 ELISA Kit (orb1212211)

Description	Mouse Glutathione Peroxidase 3, GPX3 ELISA Kit
Reactivity	Mouse
Range	1-400mU/mL
Tested Applications	ELISA
Target	GPX3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,424
Uniprot ID	P46412
Solubility (25°C)	480mU/mL
Sensitivity	0.56mU/mL
Entrez	14778
Expiration Date	Please enquire.